



synedra



information technologies

IHE Integration Statement

synedra AIM Version Prometheus

Date: 2007-05-02
Vendor: synedra information technologies GmbH

IHE Integration Statement		Date: 2007-05-02
Vendor	Product Name	Version
synedra information technologies GmbH	synedra AIM	Prometheus
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Performed Procedure Step Manager	None
Scheduled Workflow	Image Manager/Archive	Instance Availability Notification PPS Exception Management
Consistent Time	Time Server	None
Consistent Time	Time Client	None
<p>Internet address for vendor's IHE information: http://www.synedra.com/de/ihe.php</p>		
<p>Links to Standards Conformance Statements for the implementation</p>		
DICOM:	http://www.synedra.com/downloads.php	
<p>Links to general information on IHE</p>		
In North America : http://www.ihe.net	In Europe : http://www.ihe-europe.org	In Japan : http://www.jira-net.or.jp/ihe-j