



synedra




information technologies

IHE Integration Statement

synedra AIM Version 2.7 „Perseus“

Date: 2009-04-30

Vendor: synedra information technologies GmbH



Vendor	Product Name	Version
synedra information technologies GmbH	synedra AIM	2.7 „Perseus“

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
ARI-Access to Radiology Information	Image Display	None
ARI-Access to Radiology Information	Image Manager/Archive	None
ARI-Access to Radiology Information	Report Reader	None
Consistent Time	Time Client	None
Consistent Time	Time Server	None
Portable Data for Imaging	Image Display	None
Portable Data for Imaging	Portable Media Creator	None
Portable Data for Imaging	Report Reader	None
Scheduled Workflow	Image Display	None
Scheduled Workflow	Image Manager/Archive	None
Scheduled Workflow	Performed Procedure Step Manager	None
Simple Image and Numeric Report	Report Reader	None

Internet address for vendor's IHE Information

<http://www.synedra.com/de/ihe.php>

Links to Standard Conformance Statements for the implementation

DICOM: <http://www.synedra.com/downloads.php>

Links to general information on IHE

- In North America: <http://www.ihe.net>
- In Europe: <http://www.ihe-europe.org>
- In Japan: <http://www.jira-net.or.jp/ihe-j>