



synedra




information technologies

# IHE Integration Statement

synedra AIM Version 3.0 „Pegasos“

**Date:** 2011-04-18

**Vendor:** synedra information technologies GmbH



Vendor	Product Name	Version
synedra information technologies GmbH	synedra AIM	3.0 „Pegasos“

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented	Tested <sup>1</sup>
ARI-Access to Radiology Information	Image Display		2011
ARI-Access to Radiology Information	Image Manager/Archive		2009
ARI-Access to Radiology Information	Report Reader		2009
ATNA-Audit Trail and Node Authentication	Secure Application		2011
CT-Consistent Time	Time Client		2011
CT-Consistent Time	Time Server		2011
CPI-Consistent Presentation of Images	Image Display		2011
PDI-Portable Data for Imaging	Image Display		2008
PDI-Portable Data for Imaging	Portable Media Creator		2008
PDI-Portable Data for Imaging	Report Reader		2009
PDQ-Patient Demographics Query	Patient Demographics Consumer		2011
PIX-Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer		2011
SWF-Scheduled Workflow	Image Display		2009
SWF-Scheduled Workflow	Image Manager/Archive		2009
SWF-Scheduled Workflow	Performed Procedure Step Manager		2009
SINR-Simple Image and Numeric Report	Report Reader		2009
XDS.b-Cross-Enterprise Document Sharing	Document Consumer		2011
XDS.b-Cross-Enterprise Document Sharing	Integrated Document Source/Repository		2011
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Document Repository		2011
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer		2011

<sup>1</sup> Not all of our modules are tested every year, since they alter only marginally. For the sake of clarity we display here the year of the most recent test at a Connectathon for each Profile/Actor pair.

### **Internet adress for vendor´s IHE Information**

<http://www.synedra.com/de/ihe.php>

### **Links to Standard Conformance Statements for the implementation**

DICOM: <http://www.synedra.com/downloads.php>

### **Links to general information on IHE**

- In North America: [www.ihe.net](http://www.ihe.net)
- In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)
- In Japan: [www.ihe-j.org](http://www.ihe-j.org)