



synedra





information technologies

# IHE Integration Statement

synedra AIM Version 3.3 „Eos“

**Date:** 2013-04-22

**Vendor:** synedra information technologies GmbH



<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
synedra information technologies GmbH	synedra AIM	3.3 „Eos“

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Audit Trail and Node Authentication	Secure Application	
Consistent Time	Time Client	
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	Set of DICOM Instances
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	
Cross-Enterprise Document Sharing for Imaging	Document Repository	
Cross-Enterprise Document Sharing	Document Consumer	
Cross-Enterprise Document Sharing	Integrated Document Source/Repository	
Patient Demographics Query	Patient Demographics Consumer	
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	

**Internet adress for vendor´s IHE Information**

<http://www.synedra.com/de/ihe.php>

**Links to Standard Conformance Statements for the implementation**

DICOM: [http://www.synedra.com/downloads2/DICOM\\_Conformance\\_Statement\\_synedra\\_AIM.xml.pdf](http://www.synedra.com/downloads2/DICOM_Conformance_Statement_synedra_AIM.xml.pdf)

**Links to general information on IHE**

- In North America: [www.ihe.net](http://www.ihe.net)
- In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)
- In Japan: [www.ihe-j.org](http://www.ihe-j.org)