


IHE Integration Statement

synedra AIM Version 15 „Odysseus“

Date: 2015-05-20

Vendor: synedra information technologies GmbH



IHE Integration Statement

Date: 2015-05-20

Vendor	Product Name	Version
synedra information technologies GmbH	synedra AIM	15 „Odysseus“

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
ATNA-Audit Trail and Node Authentication	Audit Record Repository	None
Consistent TimelMaging Object Change Management	Time Client	None
Image Object Change Management	Image Manager/Archive	None
SWF.b-Scheduled Workflow.b	Performed Procedure Step Manager	None
SWF.b-Scheduled Workflow.b	Image Manager/Archive	None

Internet address for vendor's IHE Information

<http://www.synedra.com/de/ihe.php>

Links to Standard Conformance Statements for the implementation

DICOM: http://www.synedra.com/downloads2/DICOM_Conformance_Statement_synedra_AIM.xml.pdf

Links to general information on IHE

- In North America: www.ihe.net
- In Europe: www.ihe-europe.net
- In Japan: www.ihe-j.org