


IHE Integration Statement

synedra AIM Version 16 „Hermes“

Date: 2016-05-25

Vendor: synedra information technologies GmbH



Vendor	Product Name	Version
synedra information technologies GmbH	synedra AIM	16 „Hermes“

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
ATNA-Audit Trail and Node Authentication	Secure Application	None
ATNA-Audit Trail and Node Authentication	Audit Record Repository	None
CT-Consistent Time	Time Client	None
PDQ-Patient Demographics Query	Patient Demographics Consumer	None
PDQV3-Patient Demographics Query HL7 V3	Patient Demographics Consumer	None
PDQV3-Patient Demographics Query HL7 V3	Patient Demographics Consumer	Continuation Option
PIX-Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
PIXV3-Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
PIXV3-Patient Identifier Cross-Reference HL7 V3	Patient Identity Source	None
SWF.b-Scheduled Workflow.b	Performed Procedure Step Manager	None
SWF.b-Scheduled Workflow.b	Image Manager/Archive	None
XDR-Cross-Enterprise Document Reliable Interchange	Document Administrator	None
XDR-Cross-Enterprise Document Reliable Interchange	Document Administrator	Delete Metadata
XDR-Cross-Enterprise Document Reliable Interchange	Document Administrator	Update Metadata
XDR-Cross-Enterprise Document Reliable Interchange	Document Source	None
XDS.b-Cross-Enterprise Document Sharing	Document Administrator	Update Metadata
XDS.b-Cross-Enterprise Document Sharing	Document Administrator	None
XDS.b-Cross-Enterprise Document Sharing	Document Administrator	Delete Metadata
XDS.b-Cross-Enterprise Document Sharing	Document Consumer	None
XDS.b-Cross-Enterprise Document Sharing	Document Repository	None
XDS.b-Cross-Enterprise Document Sharing	Document Source	None
XDS.b-Cross-Enterprise Document Sharing	Integrated Document	None

Integration Profiles Implemented	Actors Implemented	Options Implemented
	Source/Repository	
XDS.b-Cross-Enterprise Document Sharing	Patient Identity Source	None
XDS.b-Cross-Enterprise Document Sharing	Patient Identity Source	Patient Identity Feed (HL7 V3)
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Document Repository	None
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	PDF Report
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	None
XDS-I.b-Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	Set of DICOM Instances
XDS-SD-Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	None
XDS-SD-Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	View Option
XDS-I.b-Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	Document Import Option
XUA-Cross Enterprise User Assertion	X-Service Provider	PurposeOfUse
XUA-Cross Enterprise User Assertion	X-Service Provider	None
XUA-Cross Enterprise User Assertion	X-Service Provider	Subject-Role
XUA-Cross Enterprise User Assertion	X-Service User	None
XUA-Cross Enterprise User Assertion	X-Service User	Subject-Role
XUA-Cross Enterprise User Assertion	X-Service User	PurposeOfUse

Internet address for vendor's IHE Information

<http://www.synedra.com/de/loesungen/ihe/>

Links to Standard Conformance Statements for the implementation

DICOM:

http://www.synedra.com/de/downloads/Dokumente_Tools/Dicom_Conformance_Statement.pdf

Links to general information on IHE

- In North America: www.ihe.net
- In Europe: www.ihe-europe.net
- In Japan: www.ihe-j.org